

Opravný lístek diplomové práce PŘF UK

Vakuolární aspartátová proteasa kvasinky *Candida albicans*

Bc. Dominik Hirko, 2024

1. Doplnění informací k fenotypovým testům

Interestingly, *C. albicans* KC281-*URA3* displayed limited growth under most conditions (Figures 39-46). This may be attributed to a technical issue with the replicator, potentially exceeding the optimal temperature for yeast viability, thus killing cells upon application to certain media. Consequently, it is imperative to repeat these experiments and re-evaluate the growth of *C. albicans* KC281-*URA3* under the specified conditions.